OPHTHALMIC INSTRUMENT TABLE BINGO 001




MediWorks®


- Easy to install(right or left)
- Including chair
- Electric Magnetic lock
- Economic
- Space-saving
- Morden design

OPHTHALMIC INSTRUMENT TABLE BINGO 001


This product is a combination of a table and chair that offers a functional use of the refraction instruments during basic examination (Refractometer , Slit Lamp, Noncontact Tonometer, Javal keratometer, Chart Projector). Ophthalmic instrument tables are designed for a faster and more comfortable use by doctor with ophthalmology instruments in examination rooms and hospitals.







SPECIFICATIONS	
Chair height	450-630mm
Loading for chair	150kg
Table size	880*400mm
Loading for table	40kg
Table height	725mm
Table lock	Magnetic
Power supply	100-230V AC/50Hz
Fuse	8A 250V
Lamp	0-240V AC/50Hz
Weight	210kg
Figure size	1660*1900*1800mm



Shanghai MediWorks Precision Instruments Co., Ltd.
Add:Building 7, Ming Pu Plaza, No. 3279, San Lu Rd, Min Hang District, Shanghai, 201100, china
Tel:+86-21-5426042154260423 Fax:+86-21-54260425
Email:marketing@mediworks.biz international@mediworks.biz