
# POLYPROPYLENE

## BLUE



| | |
|---|---|
| **Type** | Monofilament |
| **Composition** | Polypropylene, a polymer of propylene |
| **Coating** | None |
| **Colour** | Blue |
| **Tissue reaction** | Minimal |
| **Absorption** | Non absorbable |
| **Presentation** | Box with 12 sutures |

*For quality & safety*

www.sutures.be

# POLYPROPYLENE

BLUE

## CHARACTERISTICS

| | |
|---|---|
| ☑ High tensile strength, non absorbable | Excellent and permanent tissue support |
| ☑ Synthetic | Minimal tissue reaction |
| ☑ Monofilament with microscopic uniform diameter | Soft passage through the tissues without sawing, tissue drag and trauma<br>Absence of capillarity (no support of bacterial growth)<br>Excellent knot tie down |
| ☑ Pliability and softness | Good handling properties<br>Easy and secure knotting |
| ☑ Hermetically sealed packing | Guaranteed seal and product sterility |

***INDICATIONS***

General soft tissue closing and/or ligation

especially in

- Ophthalmic surgery
- Plastic surgery
- Skin closure
- Orthopaedics
- Hernias (also as a mesh)
- Gastrointestinal surgery
- General surgery
- Gynaecology
- Obstetrics

***CONTRAINDICATIONS***


- Cardiovascular surgery
- Neurosurgery
- Microsurgery


# POLYPROPYLENE


BLUE

| Side cutting spatula double armed | | Metric | USP | Length | | Order number |
|---|---|---|---|---|---|---|
| 3/8 circle  | 2 x DSP 6.4 mm 200 µm  | 0.2 | 10/0 | 30 cm | *blue* | OPP105 |
| | | 0.3 | 9/0 | 30 cm | *blue* | OPP95 |

| Side cutting spatula | | Metric | USP | Length | | Order number |
|---|---|---|---|---|---|---|
| 3/8 circle  | DSP 6 mm 200 µm  | 0.5 | 7/0 | 45 cm | *blue* | OPP71-1 |

| Side cutting spatula double armed | | Metric | USP | Length | | Order number |
|---|---|---|---|---|---|---|
| 1/4 circle  | 2 x VSP 8 mm 300 µm  | 0.7 | 6/0 | 45 cm | *blue* | OPP62 |

| Side cutting spatula | | Metric | USP | Length | | Order number |
|---|---|---|---|---|---|---|
| 1/4 circle  | VSP 16 mm 150 µm  | 0.2 | 10/0 | 10 cm | *blue* | OPP103 |

| Side cutting spatula loop | | Metric | USP | Length | | Order number |
|---|---|---|---|---|---|---|
| 1/4 circle  | VSP 15.3 mm 150 µm  | 0.2 | 10/0 | 20 cm | *blue* | OPP113L-2 |
| | VSP 16 mm 150 µm  | 0.2 | 10/0 | 20 cm | *blue* | OPP103L-2 |



www.sutures.be

# POLYPROPYLENE

BLUE

| Side cutting spatula double armed | | Metric | USP | Length | | Order number |
|---|---|---|---|---|---|---|
| Straight | 2 x GSP 16.15 mm<br>150 µm  | 0.2 | 10/0 | 20 cm | *blue* | OPP101 |
| | | 0.3 | 9/0 | 20 cm | *blue* | OPP91 |

| Side cutting spatula loop | | Metric | USP | Length | | Order number |
|---|---|---|---|---|---|---|
| Straight | GSP 16.15 mm<br>150 µm  | 0.2 | 10/0 | 15 cm | *blue* | OPP110 |

| Side cutting spatula double armed | | Metric | USP | Length | | Order number |
|---|---|---|---|---|---|---|
| Straight<br>3/8 circle  | GSP 16 mm<br>150 µm<br>DSP 4 mm<br>100 µm  | 0.2 | 10/0 | 10 cm | *blue* | OPP102 |



www.sutures.be