
SILK



	Silk Braided	*Silk Virgin*
Type	Braided multifilament	Twisted multifilament
Composition	Braided fibres from the cocoon of the silkworm	Twisted fibres from the cocoon of the silkworm
Coating	Silicone	None
Colour	Black	Blue
Tissue reaction	Moderate	
Absorption	Silk suture elicits an initial inflammatory reaction in tissues, which is followed by gradual encapsulation of the suture by fibrous connective tissues.	
Presentation	Box with 12 sutures	

For quality & safety

smi
www.sutures.be

SILK

CHARACTERISTICS	BENEFITS
☑ High tensile strength, non absorbable	Good and extended tissue support for up to 3 month
☑ Braided structure	Excellent handling properties High flexibility High tensile strength Excellent knot security
☑ Coated multifilament	Soft passage through the tissues with minimal sawing, tissue drag and trauma Good knot tie down / adjustability Reduced capillary action
☑ Hermetically sealed packing	Guaranteed seal and product sterility

INDICATIONS

General soft tissue closing and/or ligation

especially in

- General surgery
- Skin closure
- Ophthalmic surgery
- Oral surgery
- Gastrointestinal surgery
- Gynaecology
- Obstetrics

CONTRAINDICATIONS

- Cardiovascular surgery
- Neurosurgery
- Microsurgery
- Known allergies or sensitivities to silk

SILK
BRAIDED
BLACK

Side cutting spatula double armed		Metric	USP	Length		Order number
3/8 circle	2 x DSP 6.5 mm 200 µm	0.4	8/0	45 cm	*black*	OS81B
	2 x DSP 6.5 mm 240 µm	0.5	7/0	45 cm	*black*	OS71
		0.7	6/0	45 cm	*black*	OS61
	2 x DSP 8 mm 350 µm	1	5/0	30 cm	*black*	OS52

Side cutting spatula double armed		Metric	USP	Length		Order number
1/4 circle	2 x VSP 8 mm 250 µm	0.7	6/0	45 cm	*black*	OS62
	2 x VSP 10 mm 350 µm	1	5/0	45 cm	*black*	OS51

Reverse cutting		Metric	USP	Length		Order number
3/8 circle	DS 12 mm	0.7	6/0	45 cm	*black*	8071512
		1	5/0	45 cm	*black*	8101512
		1.5	4/0	75 cm	*black*	8151512
	DS 16 mm	1.5	4/0	75 cm	*black*	8151516
	DS 19 mm	1.5	4/0	75 cm	*black*	8151519

Reverse cutting double armed		Metric	USP	Length		Order number
3/8 circle	2 x DS 12 mm	0.7	6/0	45 cm	*black*	8071512-2

www.sutures.be

SILK VIRGIN BLUE

	Side cutting spatula double armed	Metric	USP	Length		Order number
3/8 circle 	2 x DSP 6.2 mm 150 µm	0.4	8/0	45 cm	*blue*	OS82
	2 x DSP 6.5 mm 200 µm	0.4	8/0	45 cm	*blue*	OS81

www.sutures.be